

U¹³C Plant Materials Expand Probing Capacities in Complex Ecosystems

New Developments in Stable Isotope Applications



Ton Gorissen¹, Ries de Visser¹, Sara N. de Visser² & Bernd P. Freyermann²

¹ IsoLife bv, P.O.Box 349, 6700 AH Wageningen, The Netherlands;

² Centre for Ecological and Evolutionary Studies, University of Groningen, P.O.Box 14, 9750 AA Haren, The Netherlands.

Uniformly (>98 at%) ¹³C-Labelled Higher Plants

Stable isotopes like ¹³C and ¹⁵N have been used in ecology for decades in tracer studies of plant materials in organic environments like natural soils or human gut systems. They are used to trace C & N or to reveal trophic levels in complex food webs.

Although De Visser *et al.* (2008) successfully identified trophic levels in natural abundance studies, low isotope-enriched plant materials have drawbacks:

1. Low resolution caused by dilution of the tracer in the natural background of organic matter pools;
2. Inadequacy for modern techniques like SIP requiring a high, uniform enrichment (e.g. U¹³C-cellulose).

De Visser *et al.* 2008. Trophic inter-actions among invertebrates in termitaria in the African savanna: a stable isotope approach. *Ecol. Entomol.* 33: 758.

Such U¹³C-labelled plant materials are only produced since 2005 by **IsoLife**. They open up new ways in studying trophic relations, decomposition of plant materials, and functional relationships in complex ecosystems (see examples A). Moreover, they improve our measurement techniques (see examples B&C).



Uniform labelling of plants with stable isotopes is like ringing birds.

Not a few, but all.

IsoLife's goal is to supply Life Science research with ¹³C- and ¹⁵N-labelled plants and plant components.

Feel free to contact us for information:

P: ++31 6 47784735 (mobile)

E: info@isolife.nl


W: www.isolife.nl

A) New Tools in Trophic Studies

i) Stable Isotope Probing (SIP) enables detection of active organisms growing on ¹³C-labelled substrates – either primary consumers or predators by using density gradient centrifugation to separate ¹²C- from ¹³C-DNA or -RNA. Subsequent sequencing of the ¹³C-bands will reveal the active, functional species in a food web. **Trophic interactions** in food webs can thus be studied with a new focus.

Example: U¹³C-labelled starch enables identifying gut bacteria. (Kovatcheva *et al.* 2009. Linking phylogenetic identities of bacteria to starch fermentation. *Environm. Microbiol.* 11: 914).

AWARD WINNING PUBLICATION.



ii) Decomposition studies with U¹³C-labelled lignin or cellulose can yield detailed information about the fraction respired, the distribution among fractions, and transformations. Moreover, the *same* experiments may yield information on active organisms by using SIP (i).

Example: Waldrop & Harden. 2008. Interactive effects of wildfire and permafrost on microbial communities and soil processes. *Global Change Biology* 14: 1.

B) New Internal Standards in MS

¹³C-Internal standards function as absolute yardsticks in MS analyses. Degradation of metabolites during the analyses and variations in sample processing are accounted for by spiking your samples with a U¹³C-labelled standard during processing.

Example: Erk *et al.* 2009. A novel method for the quantification of quinic acid in food. *Journal of AOAC International* 92: 730.

C) New Detection of Trace Gasses

Trace gas detection based on laser systems enable to monitor real time gas emission of trace-gasses like methane. Measurement of **¹³C-labelled gasses** in combination with uniform ¹³C-labelling of plants provide new powerful tools.

Example: U¹³C-labelled sage (*Salvia officinalis* L.) in a study of aerobic methane emission by plants. (Dueck *et al.* 2007. No evidence for substantial aerobic methane emission by terrestrial plants: A ¹³C-labelling approach. *New Phytol.* 175: 29).

